


Local and Segmental Therapy with PT2000



copyright by S+B medVET

S + B medVET
Systeme + Beratung

Area of application	Shoulder	Elbow	Spine	Hip	Knee	For muscular atrophy, to build muscle mass
Segmental pad position	Apply lateral to vertebrae C4 to C7 on affected side	Apply beside C6 to T1 and/or T1 to T4, either bilaterally, left and right of the vertebrae, or unilaterally.	Apply bilaterally, to the left and right of the painful areas of the spine as well as unilaterally, cranial and caudal to the affected area.	Apply beside L2 to L5 and/or T10 to L1, either bilaterally, to the left and right of the vertebrae, and/or unilaterally on the affected side. Warning: Do not place pads on bone points.	Apply to lumbar muscles from L3 caudally to S1 on the affected side	In atrophy of neural origin, the following additional treatments can be administered as needed: Front legs: C4 - T2. Quadriceps: L3, L4. Hamstrings (semitendinous and semimembranous muscles) + distal hind limb L6, S1. Variant: Apply one pad beside the corresponding vertebra and the other on the affected muscle.
Program	Begin with vibration pulse, then switch to program 2	Begin with vibration pulse, then switch to program 2	Begin with vibration pulse, then switch to program 2	Begin with vibration pulse, then switch to program 2	Begin with vibration pulse, then switch to program 2	Begin with vibration pulse, then switch to program 2
Pulse modulation	ON	ON	ON	ON	ON	ON
Local pad position	Apply medial and lateral to joint or diagonally across joint	Apply above and below, diagonally across, or on one side of the joint	Apply bilaterally, to the left and right of the painful areas of the spine as well as unilaterally, cranial and caudal to the affected area.	Apply cranial and caudal to hip joint, or diagonally across the joint, from above to below. Do not place pads on bone points.	Apply above and below, diagonally across, or on one side of the joint	Apply pad near to the origin and attachment of the affected muscle or to a functional unit (extensor or flexor - never treat both in a single session) while the pad is still on the belly of the muscle (not on the tendon or a bone point).
Program	Begin with vibration pulse, then switch to program 5 or 1	Begin with vibration pulse, then switch to program 5	Begin with vibration pulse, then switch to program 2	Begin with vibration pulse, then switch to program 5 or 3	Begin with vibration pulse, then switch to program 5	Begin briefly with vibration pulse, then switch to tapping pulse (strong but steady stimulation) or kneading pulse (slowly rising pulse - may take some getting used to). Increase intensity until contraction becomes visible.
Pulse modulation	ON	ON	ON	ON	ON	OFF!